

## Product Datasheet

### Tazemetostat hydrobromide (orb1963091)

|                        |   |
|------------------------|---|
| <b>Description</b>     | Tazemetostat hydrobromide   |
| <b>Form/Appearance</b> | Lyophilized   |
| <b>Storage</b>         | Store lyophilized at -20°C, keep desiccated. In lyophilized form, the chemical is stable for 36 months. In solution, store at -20°C and use within 3 months to prevent loss of potency. Aliquot to avoid multiple freeze/thaw cycles. |
| <b>Note</b>            | For research use only   |
| <b>Formula</b>         | $C_{34}H_{45}BrN_4O_4$  |
| <b>Purity</b>          | > 98%   |
| <b>MW</b>              | 653.65  |
| <b>CAS Number</b>      | 1467052-75-0  |
| <b>SMILES</b>          | <chem>O=C(C1=CC(C2=CC=C(CN3CCOCC3)C=C2)=CC(N(CC)C4CCOCC4)=C1C)NCC5=C(C)C=C(C)NC5=O.[H]Br</chem>   |
| <b>Expiration Date</b> | 12 months from date of receipt.   |