

Product Datasheet

Tazemetostat hydrobromide (orb1963092)

Description	Tazemetostat hydrobromide
Form/Appearance	Lyophilized
Storage	Store lyophilized at -20°C, keep desiccated. In lyophilized form, the chemical is stable for 36 months. In solution, store at -20°C and use within 3 months to prevent loss of potency. Aliquot to avoid multiple freeze/thaw cycles.
Note	For research use only
Formula	$C_{34}H_{45}BrN_4O_4$
Purity	> 98%
MW	653.65
CAS Number	1467052-75-0
SMILES	<chem>O=C(C1=CC(C2=CC=C(CN3CCOCC3)C=C2)=CC(N(CC)C4CCOCC4)=C1C)NCC5=C(C)C=C(C)NC5=O.[H]Br</chem>
Expiration Date	12 months from date of receipt.