

Product Datasheet

NPEPPS antibody (orb19718)

Description	Goat polyclonal antibody to NPEPPS
Species/Host	Goat
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	MP100 / NPEPPS
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:4000.WB: Approx 110kDa band observed in Human Brain (Cerebral Cortex) lysates (calculated MW of 98.5kDa according to NP_006301.2). Recommended concentration: 0.2-0.5 µg/ml.
Protein Sequence	KRDAESIHQYLLQRK
Clonality	Polyclonal
MW	98.5
Entrez	9520
Dilution Range	ELISA: 1:4000, WB: 0.2-0.5 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Human Brain (Ce...)