

Product Datasheet

Chicken anti Human IgG (H/L) Absorbed against Mouse, Goat and Rabbit IgG, conjugated with FITC (orb1972295)

Description	Chicken anti Human IgG (H/L) Absorbed against Mouse, Goat and Rabbit IgG, conjugated with FITC
Species/Host	Gallus
Conjugation	FITC
Tested Applications	IHC
Form/Appearance	Phosphate buffered saline, pH 7.2, 50% glycerol with 0.075% sodium azide
Concentration	See vial for concentration.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Source	Chickens were immunized with highly purified normal human IgG. After multiple immunizations, eggs were collected and the IgY fraction prepared. This fraction was immunoaffinity purified from an immunoabsorbant. Anti human IgG was further absorbed against mouse, goat and rabbit IgG and then conjugated to fluorescein isothiocyanate.
Hazard Information	This product is intended FOR RESEARCH USE ONLY, and FOR TESTS IN VITRO, not for use in diagnostic or therapeutic procedures involving humans or animals.
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)