
Product Datasheet

Chicken anti Human IgG (H/L) Absorbed against Mouse, Goat and Rabbit IgG, conjugated with FITC (orb1972295)

Description

Chicken anti Human IgG (H/L)

Rabbit IgG, conjugated with FITC

Species/Host

Gallus

Conjugation

Unconjugated

Tested

IHC

Applications**Form/Appearance**Phosphate buffered saline, pH 7.2,
50% glycerol with 0.075% sodium
azide**Concentration**

See vial for concentration.

StorageStorage: Product should be stored
at -20°C. Aliquot to avoid
freeze/thaw cycles**Note**

For research use only

Source

Chickens were immunized with highly purified normal human IgG. After multiple immunizations, eggs were collected and the IgY fraction prepared. This fraction was immunoaffinity purified from an immunoabsorbant. Anti human IgG was further absorbed against mouse, goat and rabbit IgG and then conjugated to fluorescein isothiocyanate.

Hazard Information

This product is intended FOR RESEARCH USE ONLY, and FOR TESTS IN VITRO, not for use in diagnostic or therapeutic procedures involving humans or animals.

Expiration Date

12 months from date of receipt.