

Product Datasheet

Chicken anti Rabbit IgG (H+L) (orb1972355)

Description	Chicken anti Rabbit IgG (H+L)
Species/Host	Gallus
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Form/Appearance	Phosphate buffered saline, pH 7.2 with 0.075% sodium azide
Concentration	See vial for concentration.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Source	Chickens were immunized with highly purified rabbit IgG isolated from normal pooled rabbit serum. After multiple immunizations in Freund's adjuvant, eggs were collected and the IgY fraction prepared. The antibodies were immunoaffinity purified from a solid phase immunoabsorbent.
Hazard Information	This product is intended FOR RESEARCH USE ONLY, and FOR TESTS IN VITRO, not for use in diagnostic or therapeutic procedures involving humans or animals.
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)