

---

## Product Datasheet

**Chicken anti Rabbit IgG (H+L), absorbed with human and mouse IgG (orb1972356)**

**Description**

Chicken anti Rabbit IgG (H+L),

IgG

**Species/Host**

Gallus

**Conjugation**

Unconjugated

**Tested**

ELISA, WB

**Applications****Form/Appearance**Phosphate buffered saline, pH 7.2  
with 0.075% sodium azide**Concentration**

See vial for concentration.

**Storage**Storage: Product should be  
stored at -20°C. Aliquot to avoid  
freeze/thaw cycles**Note**

For research use only

**Source**

Chickens were immunized with highly purified rabbit IgG isolated from normal pooled rabbit serum, in Freund's adjuvant. Eggs were collected, the IgY fraction prepared and the anti rabbit IgG immunoaffinity purified from a solid phase immunoabsorbant. Anti Rabbit IgG was further absorbed against mouse and human IgG.

**Hazard Information**

This product is intended FOR RESEARCH USE ONLY, and FOR TESTS IN VITRO, not for use in diagnostic or therapeutic procedures involving humans or animals.

**Expiration Date**

12 months from date of receipt.