

Product Datasheet

Chicken anti Rabbit IgG (H+L), conjugated with Agarose (orb1972359)

Description	Chicken anti Rabbit IgG (H+L), conjugated with Agarose
Species/Host	Gallus
Conjugation	Agarose
Form/Appearance	Phosphate buffered saline, pH 7.2 with 0.075% sodium azide
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Source	Chickens were immunized with highly purified rabbit IgG isolated from normal pooled rabbit serum. After multiple immunizations in Freund's adjuvant, eggs were collected and the IgY fraction prepared. The antibodies were immunoaffinity purified from a solid phase immunoabsorbent and conjugated to conjugated to N- hydroxysuccinimide (NHS)-activated agarose beads. Beads were washed extensively after the blocking of the residual NHS sites.
Hazard Information	This product is intended FOR RESEARCH USE ONLY, and FOR TESTS IN VITRO, not for use in diagnostic or therapeutic procedures involving humans or animals.
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)