

Product Datasheet

Chicken anti Sheep IgG (H+L) (orb1972361)

Description	Chicken anti Sheep IgG (H+L)
Species/Host	Gallus
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Form/Appearance	Phosphate buffered saline, pH 7.2 with 0.075% sodium azide
Concentration	See vial for concentration.
Storage	Storage: Product should be stored at -20°C. Aliquot to avoid freeze/thaw cycles
Note	For research use only
Source	Chickens were immunized with highly purified sheep IgG isolated from normal pooled sheep serum. After multiple immunizations in Freund's adjuvant, eggs were collected and the IgY fraction prepared. The antibodies were immunoaffinity purified from a solid phase immunoabsorbent.
Hazard Information	This product is intended FOR RESEARCH USE ONLY, and FOR TESTS IN VITRO, not for use in diagnostic or therapeutic procedures involving humans or animals.
Expiration Date	12 months from date of receipt.