

## Product Datasheet

# Chicken anti Sheep IgG (H+L), conjugated with FITC (orb1972363)

<b>Description</b>	Chicken anti Sheep IgG (H+L), conjugated with FITC
<b>Species/Host</b>	Gallus
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	IHC
<b>Form/Appearance</b>	Phosphate Buffered Saline, 50% Glycerol with 0.075% sodium azide
<b>Concentration</b>	See vial for concentration.
<b>Storage</b>	Storage: Product should be stored at -20°C. Aliquot to avoid freeze/thaw cycles
<b>Note</b>	For research use only
<b>Source</b>	Chickens were immunized with highly purified sheep IgG isolated from pooled normal sheep serum. After multiple immunizations in Freund's adjuvant, eggs were collected and the IgY fraction prepared. The antibodies were immunoaffinity purified from a solid phase immunoabsorbant and conjugated to fluorescein isothiocyanate.
<b>Hazard Information</b>	This product is intended FOR RESEARCH USE ONLY, and FOR TESTS IN VITRO, not for use in diagnostic or therapeutic procedures involving humans or animals.
<b>Expiration Date</b>	12 months from date of receipt.