

Product Datasheet

Rabbit anti Bisphenol A (orb1972503)

Description	Rabbit anti Bisphenol A
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, IAC
Immunogen	BSA-Bisphenol Valeric Acid (BVA)
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	Ig
Clonality	Polyclonal
Source	Rabbit immune serum, purified by caprylic acid extraction.
Expiration Date	12 months from date of receipt.