

## Product Datasheet

### Rabbit anti Bisphenol A (orb1972504)

<b>Description</b>	Rabbit anti Bisphenol A
<b>Species/Host</b>	Rabbit
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA, IAC
<b>Immunogen</b>	BSA-Bisphenol Valeric Acids (BVA)
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Isotype</b>	Ig
<b>Clonality</b>	Polyclonal
<b>Source</b>	Rabbit immune serum
<b>Expiration Date</b>	12 months from date of receipt.