

Product Datasheet

Anti-POGK Antibody (orb1972558)

Description	Biorbyt Anti-POGK Antibody catalog # orb1972558. Tested in WB, ICC/IF, FCM(Fixed), ELISA application...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC, ICC, WB
Immunogen	E.coli-derived human POGK recombinant protein (Position: R197-R437). Human POGK shares 97.9% amino acid (aa) sequence identity with mouse POGK.
Form/Appearance	Lyophilized
Storage	At -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freezing and thawing.
Note	For research use only
Application notes	Western blot, 0.25-0.5 µg/ml, Human, Mouse, Rat Immunocytochemistry/Immunofluorescence, 5 µg/ml, Human Flow Cytometry (Fixed), 1-3 µg/1x10 ⁶ cells, Human ELISA, 0.1-0.5 µg/ml, Human. Tested Species: In-house tested species with positive results. Other applications have not been tested. Optimal dilutions should be determined by end users. Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml
Clonality	Polyclonal
Uniprot ID	Q9P215
Expiration Date	12 months from date of receipt.