

Product Datasheet

Human CCL15 Recombinant Protein (orb1972728)

Description	The Human CCL15 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The...
Target	CCL15
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	QFTNDAETEL MMSKLPLENP VVLNSFHFAA DCCTSYISQS IPCSLMKSYF ETSSECSKPG VIFLTKKGRQ VCAKPSGPGV QDCMKKLPY SI (92)
Purity	98%
Source	Yeast
MW	10.2 kDa
Entrez	601393
Expiration Date	6 months from date of receipt.