

Product Datasheet

Feline PD-L1 Recombinant Protein (orb1972739)

Description	The Feline PD-L1 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. Th...
Target	PD-L1
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	FTITVSKDLY VVEYGSNVTM ECRFPVEEQL DLVSLIVYWE MEDKKIIQFV QGKEDLKVQH RSYSQRAQLL KDQLFLGKAA LQITNVTLED AGVYCCLIGY GGADYKRITL KVHAPYRKIN QRISVDPVTS EHELMCQAEG YPTAEVIWTN SAHQVLNGKT IISVSNMETK LFNVTSTLRI NTTANEIFYC TFLQRSSPEG NSTAELVIPE PFLVPANER (219)
Purity	98%
Source	Yeast
MW	24.6 kDa
Entrez	100127110
Expiration Date	6 months from date of receipt.