

Product Datasheet

Human CCL16 Recombinant Protein (orb1972741)

Description	The Human CCL16 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The...
Target	CCL16
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	QPKVPEWVNT PSTCCLKYYE KVLPRRLVVG YRKALNCHLP AIIFVTKRNR EVCTNPNDWW VQEYIKDPNL PLLPTRNLST VKIITAKNGQ PQLLSQ (97)
Purity	98%
Source	Yeast
MW	11.2 kDa
Entrez	601394
Expiration Date	6 months from date of receipt.