

Product Datasheet

Human CCL1 Recombinant Protein (orb1972745)

Description	The Human CCL1 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The ...
Target	CCL2
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	KSMQVPFSRC CFSFAEQEIP LRAILCYRNT SSICSNEGLI FKLKRGKEAC ALDTVGWVQR HRKMLRHCPK KKK (73)
Purity	98%
Source	Yeast
MW	8.5 kDa
Entrez	182281
Expiration Date	6 months from date of receipt.