

Product Datasheet

Human BAFF Recombinant Protein (orb1972747)

Description	The Human BAFF yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The ...
Target	BAFF
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	AVQGPEETVT QDCLQLIADS ETPTIQKGSY TFVPWLLSFK RGSALEEKEN KILVKETGYF FIYGQVLYTD KTYAMGHLIQ RKKVHVFGDE LSLVTLFRCI QNMPETLPNN SCYSAGIAKL EEGDELQLAI PRENAQISLD GDVTFFGALK LL (152)
Purity	98%
Source	Yeast
MW	17.0 kDa
Entrez	10673
Expiration Date	6 months from date of receipt.