

Product Datasheet

Human CXCL11 Recombinant Protein (orb1972753)

Description	The Human CXCL11 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. Th...
Target	CXCL11
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	FPMFKRGRCL CIGPGVKAVK VADIEKASIM YPSNNCDKIE VIITLKENKG QRCLNPKSKQ ARLIKKVER KNF (73)
Purity	98%
Source	Yeast
MW	8.3 kDa
Entrez	604852
Expiration Date	6 months from date of receipt.