

Product Datasheet

Mouse CTLA-4 Recombinant Protein (orb1972759)

Description	The Mouse CTLA-4 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. Th...
Target	CTLA-4
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	EAIQVTQPSV VCLASSHGVAS FPCEYSPSHN TDEVVRTVLR QTNDQMTEVC ATTFTEKNTV GFLDYPFCSG TFNESRVNLT IQGLRAVDTG LYLCKVELMY PPPYFVGMGN GTQIYVIDPE PCPDS (126)
Purity	98%
Source	Yeast
MW	13.9 kDa
Entrez	12477
Expiration Date	6 months from date of receipt.