

Product Datasheet

Cynomolgus Monkey IL-3 Recombinant Protein (orb1972761)

Description	The Cynomolgus Monkey IL-3 yeast-derived recombinant protein is not tagged and is naturally endotoxi...
Target	IL-3
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	APMTQTTSLK TSWAKCSNMI DEIITHLNQP PLPSPDFNNL NEEDQTILVE KNLRRSNLEA FSKAVKSLQN ASAIESILKN LPPCLPMATA APTRPPIRIT NGDRNDFRRK LKFYKLTLEN EQAQ (124)
Purity	98%
Source	Yeast
MW	14.0 kDa
Entrez	102118362
Expiration Date	6 months from date of receipt.