

Product Datasheet

Dolphin IL-24 Recombinant Protein (orb1972771)

Description	The Dolphin IL-24 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. T...
Target	IL-24
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	QEFQFGPCQV EGVVLQELWE AFQAMKDIVQ AQDNITSVQL LRKEVLQNVS EAESCYLIHA LLEFYLN TVF KNYHDKAVEF RILKSFSTLA NNFIVIVSKL QPSQEKEMFS ISESARRRFL LFQRAFKQLD REAAVTKAFG EVDILLTWME KFYQI (155)
Purity	98%
Source	Yeast
MW	18.1 kDa
Entrez	101327412
Expiration Date	6 months from date of receipt.