

Product Datasheet

Human Defensin B1 Recombinant Protein (orb1972772)

Description	The Human Defensin B1 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free...
Target	Defensin B1
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	DHYNCVSSGG QCLYSACPIF TKIQGTCYRG KAKCCK (36)
Purity	98%
Source	Yeast
MW	3.9 kDa
Entrez	1672
Expiration Date	6 months from date of receipt.