

Product Datasheet

Canine IL-16 Recombinant Protein (orb1972775)

Description	The Canine IL-16 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. Th...
Target	IL-16
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	PDLNSSTDST ASASVASDVS VDPTEATVHT VTLEKTSAGL GFSLEGGKGS LLGDKPLTIN RIFKGAASEQ SETIQPGDEI LHLAGTAMQG LTRFEAWNII KTLPDGPIMM VIRRKGLQSK GITAAGDA (128)
Purity	98%
Source	Yeast
MW	13.2 kDa
Entrez	488762
Expiration Date	6 months from date of receipt.