
Product Datasheet

Bovine IL-12 p70 Recombinant Protein (orb1972789)

Description

The Bovine IL-12 p70 yeast-derived

is naturally endotoxin-free...

Target

IL-12 p70

Form/Appearance

Lyophilized

Storage

-20°C

Note

For research use only

Protein Sequence

MWELEKNVYV VELDWYPDAP
GETVVLTCDT PEEDGITWTS
DQSSEVLGSG KTLTIQVKEF
GDAGQYTCHK GGEALSRLL
LLHKKEDGIW STDILKDQKE
PKAKSFLKCE AKDYSGHFTC
WWLTAISTDL KFSVKSSRGS
SDPRGVTCGA ALLSAEKVSL
EHREYNKYTV ECQEGSACPA
AEESLLIEVV VEAHKLKYE NYTSSFFIRD
IIKPDPPKNL QLRPLKNSRQ
VEVSWEYPDT WSTPHSYFSL
TFCVQVQGKN KREKLFMDQ
TSAKVTCHKD ANVRVQARDR
YYSSFWSEWA SVSCS (305) **LINKER**
RSLPTTTASP GRSCLDYSQN
LLRAVSNTLQ KARQTLEFYS
CTSEEIDHED ITKDKTSTVE ACLPLELATN
ESCLASRETS FITNGHCLAS
GKTSFMTTLC LRSIYEDLKM
YHVEFQAMNA KLLMDPKRQI
FLDQNMLAAI AELMQALNFD
SETVPQKPSL KELDFYKTKV KLCILLHAFR
IRAVTIDRMM SYLSSS (196)

Purity

98%

Source

Yeast

MW

57.6 kDa

Entrez[35](#), [281856](#), [281857](#), [70](#)**Expiration Date**

6 months from date of receipt.