

Product Datasheet

Human CCL14 Recombinant Protein (orb1972791)

Description	The Human CCL14 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The...
Target	CCL14
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	TKTESSSRGP YHPSECCFTY TTYKIPRQRI MDYYETNSQC SKPGIVFITK RGHSVCTNPS DKWVQDYIKD MKEN (74)
Purity	98%
Source	Yeast
MW	8.7 kDa
Entrez	6358
Expiration Date	6 months from date of receipt.