

Product Datasheet

Multi-species IL-11 (Canine, Feline, Ferret, Mink) Recombinant Protein (orb1972802)

| | |
|-------------------------|--|
| Description | The Multi-species IL-11 (Canine, Feline, Ferret, Mink) yeast-derived recombinant protein is not tagg... |
| Target | IL-11 |
| Form/Appearance | Lyophilized |
| Storage | -20°C |
| Note | For research use only |
| Protein Sequence | APGPPPGPPR ASPDPRAELD SAVLLTRSLL ADTRQLAAQL RDKFPADGDH NLDLPTLAM SAGALGALQL PGVLTRLRAD LFSYLRHVQW LRRAGGPSLR TLEPELGALQ ARLDRLLRRL QLLMSRLALP QAPPDPPAPP LAPPASAWGG IRAAHAILGG LHLTDWAVR GLLLLKTRL (179) |
| Purity | 98% |
| Source | Yeast |
| MW | 19.1 kDa |
| Entrez | 101085086 |
| Expiration Date | 6 months from date of receipt. |