

## Product Datasheet

# Multi-species IL-11 (Canine, Feline, Ferret, Mink) Recombinant Protein (orb1972803)

<b>Description</b>	The Multi-species IL-11 (Canine, Feline, Ferret, Mink) yeast-derived recombinant protein is not tagg...
<b>Target</b>	IL-11
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	-20°C
<b>Note</b>	For research use only
<b>Protein Sequence</b>	APGPPPGPPR ASPDPRAELD SAVLLTRSLL ADTRQLAAQL RDKFPADGDH NLDLPTLAM SAGALGALQL PGVLTRLRAD LFSYLRHVQW LRRAGGPSLR TLEPELGALQ ARLDRLLRRL QLLMSRLALP QAPPDPPAPP LAPPASAWGG IRAAHAILGG LHLTDWAVR GLLLLKTRL (179)
<b>Purity</b>	98%
<b>Source</b>	Yeast
<b>MW</b>	19.1 kDa
<b>Entrez</b>	<a href="#">101085086</a>
<b>Expiration Date</b>	6 months from date of receipt.