

Product Datasheet

Multi-species IL-11 (Canine, Feline, Ferret, Mink) Recombinant Protein (orb1972804)

Description	The Multi-species IL-11 (Canine, Feline, Ferret, Mink) yeast-derived recombinant protein is not tagg...
Target	IL-11
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	APGPPPGPPR ASPDPRAELD SAVLLTRSLL ADTRQLAAQL RDKFPADGDH NLDLPTLAM SAGALGALQL PGVLTRLRAD LFSYLRHVQW LRRAGGPSLR TLEPELGALQ ARLDRLLRRL QLLMSRLALP QAPPDPPAPP LAPPASAWGG IRAAHAILGG LHLTDWAVR GLLLLKTRL (179)
Purity	98%
Source	Yeast
MW	19.1 kDa
Entrez	101085086
Expiration Date	6 months from date of receipt.