

Product Datasheet

Mouse IL-2 Recombinant Protein (orb1972811)

Description	The Mouse IL-2 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The ...
Target	IL-2
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	APTSSSTSSS TAEAQQQQQ QQQQQHLEQ LLMDLQELLS RMENYRNLKL PRMLTFKFYL PKQATELKDL QCLEDELGPL RHVLDLTQSK SFQLEDAENF ISNIRVTVVK LKGS DNTEFC QFDDESATVV DFLRRWIAFC QSIISTSPQ (149)
Purity	98%
Source	Yeast
MW	17.2 kDa
Entrez	16183
Expiration Date	6 months from date of receipt.