

Product Datasheet

Guinea Pig IL-4 Recombinant Protein (orb1972814)

Description	The Guinea Pig IL-4 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free....
Target	IL-4
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	CNHHTLQEII QHLNTLSREK SPCAELLVTD VFADPQGPAS GDLCTAATVL HHTAYLRGPQ SCPNREGDPL YPSVLRQVFR NLRSMASQNC PVSELRQTTL KDFLENLKRI MQKRYSKCRR (120)
Purity	98%
Source	Yeast
MW	13.7 kDa
Entrez	100720403
Expiration Date	6 months from date of receipt.