

Product Datasheet

Human IL-11 Recombinant Protein (orb1972816)

Description	The Human IL-11 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The...
Target	IL-11
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	RVSPDPRAEL DSTVLLTRSL LADTRQLAAQ LRDKFPADGD HNLDSLPTLA MSAGALGALQ LPGVLTRLRA DLLSYLRHVQ WLRRAGGSSL KTLPELGTL QARLDRLRR LQLLMSRLAL PQPPDPPAP PLAPPSSAWG GIRAAHAILG GLHLTLDWAV RGLLLLKTRL (170)
Purity	98%
Source	Yeast
MW	18.4 kDa
Entrez	3589
Expiration Date	6 months from date of receipt.