

Product Datasheet

Chicken ANG-4 Recombinant Protein (orb1972825)

Description	The Chicken ANG-4 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. T...
Target	ANG-4
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	QNEGYEKFLR QHYDAKPKGR DDRYCESMMK ERKLTSPCKD VNTFIHGTTK NIRAIKGKKG SPYGENFRIS NSPFQITTCT HSGASPRPPC GYRAFKDFRY IVIACEDGWP VHFDEFISF (120)
Purity	98%
Source	Yeast
MW	13.8 kDa
Entrez	100302572
Expiration Date	6 months from date of receipt.