

Product Datasheet

Multi-species IL-5 (Bovine, Caprine, Ovine) Recombinant Protein (orb1972827)

| | |
|-------------------------|---|
| Description | The Bovine/Caprine/Ovine IL-5 yeast-derived recombinant protein is not tagged and is naturally endot... |
| Target | IL-5 |
| Form/Appearance | Lyophilized |
| Storage | -20°C |
| Note | For research use only |
| Protein Sequence | VESTMNRLVA ETLTLLSSHR TLLIGDGNLM IPTPQHTNHQ LCIEEVFQGI DTLKNQTAQG DAVKKIFQNL SLIKEYIDLQ KRKCGGERWR VKQFLDYLQV FLGVINTEWT MES (113) |
| Purity | 98% |
| Source | Yeast |
| MW | 13.0 kDa |
| Entrez | 280825 |
| Expiration Date | 6 months from date of receipt. |