

Product Datasheet

Ferret/Mink LIF Recombinant Protein (orb1972830)

Description	The Ferret/Mink Leukemia Inhibitory Factor yeast-derived recombinant protein is not tagged and is na...
Target	LIF
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	TPLPITPVNA TCATRHPCHS NLMNQIRNQL AQLNGSANAL FILYYTAQGE PFPNNLDKLC GPNVTDFFPF HANGTEKTRL VELYRIIAYL GASLGNITRD QKVLNPNALS LHSKLNATAD ILRGLLSNVL CRLCNKYHVA HVDVAYGPD SGKDVFQKKK LGCQLLGKYK QVIAVVAQAF (180)
Purity	98%
Source	Yeast
MW	19.7 kDa
Entrez	101674557
Expiration Date	6 months from date of receipt.