

## Product Datasheet

### Ferret/Mink LIF Recombinant Protein (orb1972831)

<b>Description</b>	The Ferret/Mink Leukemia Inhibitory Factor yeast-derived recombinant protein is not tagged and is na...
<b>Target</b>	LIF
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	-20°C
<b>Note</b>	For research use only
<b>Protein Sequence</b>	TPLPITPVNA TCATRHPCHS NLMNQIRNQL AQLNGSANAL FILYYTAQGE PFPNNLDKLC GPNVTDFFPF HANGTEKTRL VELYRIIAYL GASLGNITRD QKVLNPNALS LHSKLNATAD ILRGLLSNVL CRLCNKYHVA HVDVAYGPD SGKDVFQKKK LGCQLLGKYK QVIAVVAQAF (180)
<b>Purity</b>	98%
<b>Source</b>	Yeast
<b>MW</b>	19.7 kDa
<b>Entrez</b>	<a href="#">101674557</a>
<b>Expiration Date</b>	6 months from date of receipt.