

Product Datasheet

Ovine LIF Recombinant Protein (orb1972835)

Description	The Ovine LIF yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The p...
Target	LIF
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	SPLPITPVNA TCATRHPCPS NLMNQIRNQL GQLNSSANSL FILYYTAQGE PFPNNLDKLC SPNVTFDFPPF HANGTEKARL VELYRIIAYL GASLGNITRD QKVLNPHYAHS LHSKLNATAD VLRGLLSNVL CRLCSKYHVS HVDVTYGPDT SGKDIFQKKK LGCQLLGKYK QVIAVLAQAF (180)
Purity	98%
Source	Yeast
MW	19.8 kDa
Entrez	101110180
Expiration Date	6 months from date of receipt.