

Product Datasheet

Hamster CXCL10 Recombinant Protein (orb1972845)

Description	The Hamster CXCL10 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. ...
Target	CXCL10
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	IPLSRTVRCS CIKIDDRPVK PRALGKLEII PASQSCPVE IIVTMKKTEE KRCLNPESEA IKSLLKAVSQ RRSKRAS (77)
Purity	98%
Source	Yeast
MW	8.7 kDa
Entrez	101840942
Expiration Date	6 months from date of receipt.