

Product Datasheet

Feline IL-6 Recombinant Protein (orb1972860)

Description	The Feline IL-6 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The...
Target	IL-6
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	TPGPLGGDAT SNRLPLTSAD KMEELIKYIL GKISALKKEM CDNYNKCEDS KEALAENNLN LPKLAEKDGC FQSGFNQETC LTRITTGLQE FQIYKFLQD KYEGDKENAK SVYTSTNVLL QMLKRKGKNQ DEVTIPVPTV EVGLQAKLQS QEEWLRHTTI HLTLRRLEDF LQFSLRAVRI M (181)
Purity	98%
Source	Yeast
MW	20.6 kDa
Entrez	493687
Expiration Date	6 months from date of receipt.