

Product Datasheet

Dolphin IL-17A Recombinant Protein (orb1972866)

Description	The Dolphin IL-17A yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. ...
Target	IL-17A
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	GIIPQSPGC PHVEDKNFPQ YVRLNLNILN RNVNSRRPSD YYKRSTSPWT LHRNEDSERY PSVIWEAKCS HLGCVNAEGK VDHHMNSVPI QQEILVLRRE PQHCPHSFRL EKMLVAVGCT CVTPIVHHMA (130)
Purity	98%
Source	Yeast
MW	14.9 kDa
Entrez	101322970
Expiration Date	6 months from date of receipt.