

Product Datasheet

Rabbit IL-1 beta Recombinant Protein (orb1972868)

Description	The Rabbit IL-1 beta yeast-derived recombinant protein is not tagged and is naturally endotoxin-free...
Target	IL-1 beta
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	AVRSLHCRLQ DAQQKSLVLS GTYELKALHL NAENLNQQVV FSMSFVQGEE SNDKIPVALG LRGKNLYLSC VMKDDKPTLQ LESVDPNRYK KKKMEKRFVF NKIEIKDKLE FESAQFPNWY ISTSQTEYMP VFLGNNSGGQ DLIDFSMEFV SS (152)
Purity	98%
Source	Yeast
MW	17.4 kDa
Entrez	100008990
Expiration Date	6 months from date of receipt.