

Product Datasheet

Equine LIF Recombinant Protein (orb1972871)

Description	The Equine Leukemia Inhibitory Factor yeast-derived recombinant protein is not tagged and is natural...
Target	LIF
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	SPLPITPDNA TCATRHPCHS NLMDQIRNQL VQLNSSANAL FILYYTAQGE PFPSNLDKLC RPDVTFDFPPF HANGKEKARL VELYRIIAYL GASLGNITRD QKVLNPNALS LHSKLNATSD TMRGLLSNVL CHLCSKYHVA NVDVAYGPDV SNKDVVFQKKK LGCQLLGKYK QVIAVVAQAF QTGDLEAYSD PRS (193)
Purity	98%
Source	Yeast
MW	21.3 kDa
Entrez	100629131
Expiration Date	6 months from date of receipt.