

Product Datasheet

Human CCL13 Recombinant Protein (orb1972874)

Description	The Human CCL13 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The...
Target	CCL14
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	QPDALNVPST CCFTFSSKKI SLQRLKSYVI TTSRCPQKAV IFRTKLGKEI CADPKEKWVQ NYMKHLGRKA HTLKT (75)
Purity	98%
Source	Yeast
MW	8.6 kDa
Entrez	6357
Expiration Date	6 months from date of receipt.