

Product Datasheet

Bovine TIM-3 Recombinant Protein (orb1972885)

Description	The Bovine TIM-3 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. Th...
Target	TIM-3
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	AYVSQVGQNA DLPCTYSPAT TENLVPVCWG KGPCPVFECY SLVLRDGRN VTYQTSSRYL LKRDLHKGDV TLTIKNVTLA DSGTYCCRIQ FPGLMNRKS NLELIIPAK VTPAWTPWRD ITTAAPRMLT TKGPVSETRT LKTLHDKNQT EISTLATELQ DMGATTRTG (169)
Purity	98%
Source	Yeast
MW	18.8 kDa
Entrez	768047
Expiration Date	6 months from date of receipt.