

Product Datasheet

Hamster CCL28 Recombinant Protein (orb1972893)

Description	The Hamster CCL28 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. T...
Target	CCL28
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	ILPIASSCCT EVSHHISRRL LERVNSCSIQ RADGDCDLAA VILHVKRRRI CVSPHNQTLK HWMRVSEEKK NGKENPCHKK LHSRRDRKSH IQRNHGKQGH NAPR (104)
Purity	98%
Source	Yeast
MW	12.0 kDa
Entrez	101843973
Expiration Date	6 months from date of receipt.