

Product Datasheet

Guinea Pig CCL20 Recombinant Protein (orb1972894)

Description	The Guinea Pig CCL20 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free...
Target	CCL20
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	ASSFDCCLRY TRHVLSPKLI VGFTQQLASE ACDIDAIIFH TKKRLSVCAD PPKGWVKRAV LVLSQKVKKM (70)
Purity	98%
Source	Yeast
MW	7.9 kDa
Entrez	100714505
Expiration Date	6 months from date of receipt.