

Product Datasheet

Canine CTLA-4 Recombinant Protein (orb1972904)

Description	The Canine CTLA-4 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. T...
Target	CTLA-4
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	KGMHVAQPAV VCLASSRGVAS FVCEYGSSGN AAEVRVTVLR QAGSQMTEVC AATYTVEDEL AFLDDSTCTG TSSGNKVNLT IQGLRAMDTG LYICKVELMY PPPPYYVGMGN GTQIYVIDPE PCPDSD (126)
Purity	98%
Source	Yeast
MW	13.4 kDa
Entrez	403696
Expiration Date	6 months from date of receipt.