

Product Datasheet

Human CCL8 Recombinant Protein (orb1972912)

Description	The Human CCL8 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The ...
Target	CCL8
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	QPDSVSIPIT CCFNVINRKI PIQRLESYTR ITNIQCPKEA VIFKTKRGKE VCADPKERWV RDSMKHLDQI FQNLKP (76)
Purity	98%
Source	Yeast
MW	8.9 kDa
Entrez	6355
Expiration Date	6 months from date of receipt.