

Product Datasheet

Bee Apidaecin 1A Recombinant Protein (orb1972920)

| | |
|-------------------------|---|
| Description | The Bee Apidaecin 1A yeast-derived recombinant protein is not tagged and is naturally endotoxin-free... |
| Target | Apidaecin |
| Form/Appearance | Lyophilized |
| Storage | -20°C |
| Note | For research use only |
| Protein Sequence | GNNRPVYIPQ PRPPHPR (17) |
| Purity | 98% |
| Source | Yeast |
| MW | 2.0 kDa |
| Entrez | 406140 |
| Expiration Date | 6 months from date of receipt. |