

Product Datasheet

Hamster IL-13 Recombinant Protein (orb1972928)

Description	The Hamster IL-13 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. T...
Target	IL-13
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	APGPVPRSMS PPVALRELIE ELSNITQDQR APLCNGSMVW SVDLTNGGFC AALDSLTNIS SCRPIYKTQR ILNGLCVRRA SAGVSSLPDT KVEVVQFIVT LLSYSRQLYR HGKFP (115)
Purity	98%
Source	Yeast
MW	12.6 kDa
Entrez	101843388
Expiration Date	6 months from date of receipt.